

## TN-1358

# Determination of Semaglutide and Tirzepatide in Plasma

Jiao Wang<sup>1</sup>, Liuliu Wu<sup>1</sup> and Bryan Tackett, PhD<sup>2</sup>

<sup>1</sup>Tianjin Bona Agela Co., Ltd, 179 South Street, West Zone, Tianjin Economic-Technological Development Area (TEDA), Tianjin, China, 300462

<sup>2</sup>Phenomenex Inc., 411 Madrid Ave., Torrance, CA 90501 USA



### Introduction

Glucagon-like peptide-1 (GLP-1) agonists are considered one of the most promising classes of diabetes drugs on the market. Semaglutide and Tirzepatide, as next-generation GLP-1 drugs, not only stimulate insulin release in the body effectively controlling blood sugar levels, but also inhibit gastrointestinal motility and increase satiety. This technical note establishes a method for determining Semaglutide and Tirzepatide drugs in plasma, providing a scientific basis for further research on drug safety and therapeutic efficacy.

### Sample Preparation

Step	Description
<b>Sample Pre-treatment:</b>	Combine 200 µL of EDTA anticoagulated bovine plasma and 400 µL of Methanol for protein precipitation. Centrifuge at 15,000 g for 5 minutes to obtain 500 µL supernatant. Add 400 µL of Water, then vortex and use SPE to further clean up the sample.
<b>Condition:</b>	Peptide-3-MW, 5 mg/1 mL 96-well plate (Agela) with 200 µL of Methanol, then 200 µL Water.
<b>Load:</b>	Pre-treated samples into wells.
<b>Wash:</b>	200 µL Water.
<b>Elute:</b>	50 µL elution solvent (5 % Formic Acid in Ethanol / Water (4:1, v/v)) twice. Vortex to mix.
<b>Inject:</b>	5 µL

**Table 1.** Semaglutide and Tirzepatide Structural Information.

Name	Molecular Weight (Da)	Structural Formula
Semaglutide	4113.57	C <sub>187</sub> H <sub>291</sub> N <sub>45</sub> O <sub>59</sub>
Tirzepatide	4813.45	C <sub>225</sub> H <sub>348</sub> N <sub>48</sub> O <sub>68</sub>

### LC Conditions

**Column:** Aeris™ 2.6 µm Peptide XB-C18  
**Dimensions:** 100 x 2.1 mm  
**Part No.:** [00D-4505-AN](#)  
**Mobile Phase:** A: 0.1 % Formic Acid  
 B: 0.1 % Formic Acid in Acetonitrile  
**Gradient:**

Time (min)	%B
0	30
0.5	30
3	65
3.5	65
4	98
5.5	98
5.6	30
7	30

**Flow Rate:** 0.3 mL/min  
**Injection Volume:** 5 µL  
**Temperature:** 40 °C  
**LC System:** Shimadzu® LC-20AD  
**Detection:** MRM  
**Detector:** SCIEX® 6500 Triple Quad™

### MS Conditions

**Ion Source:** ESI  
**Scan Mode:** MRM Positive and Negative Mode  
**Source Temperature:** 450 °C  
**GS1:** 60 psi  
**GS2:** 60 psi  
**CUR:** 30 psi  
**CAD:** High  
**IS:** 5500 V

**Table 2.** MRM Transitions and Parameters.

Name	Q1 (m/z)	Q3 (m/z)	DP (V)	CE (V)
Semaglutide	1029.3	1238.5	40	41
	1029.3	1110.3	40	39
	1029.3	690.2	40	39
	1029.3	960.5	40	51
Tirzepatide	1204.2	396.3	67	36
	1204.2	910.0	67	33
	1204.2	795.8	67	35




**Results and Discussion**


As shown in **Figure 1**, the linear range of the method for both compounds was 0.5 ng/mL to 200 ng/mL, with  $R^2 > 0.995$ . The limit of quantitation for both compounds was achievable at 0.5 ng/mL plasma concentration. **Figures 2 and 3** demonstrate well-defined peaks for both compounds with minimal residual carryover (**Figure 4**). Upon reaching the highest point on the curve at 200 ng/mL, the residual levels in the blank solution are below one-thousandth for Semaglutide and within five-ten-thousandths for Tirzepatide.

The spiked recovery and matrix effect were evaluated for plasma concentrations of 10 ng/mL and 100 ng/mL, respectively. The results are shown in **Table 4**. Both compounds exhibited absolute recoveries of over 80 %, however, they demonstrated some degree of matrix suppression effect.

**Figure 1.** Standard Curve for Semaglutide (Pink) and Tirzepatide (Blue).



**Figure 2.** LQD – 0.5 ng/mL Matrix Standard Quantification Ion Chromatogram for Semaglutide (Left) and Tirzepatide (Right).



**Figure 3.** Chromatogram of Matrix Standard at the Highest Point of the Curve (200 ng/mL).




Figure 4. Chromatogram of Blank Residuals after High Concentration Points on the Curve.




Table 3. Results of Spiked Recovery and Matrix Effect.

Compound	Semaglutide		Tirzepatide		
	Standard Concentration	10 ng/mL	100 ng/mL	10 ng/mL	100 ng/mL
% Recovery		82.4	95.8	81.3	93.5
% Matrix Effect		76.6	87.0	60.6	58.8

**Conclusions**

This experiment established a quantitative method for Semaglutide and Tirzepatide in plasma. The limit of quantitation (LQD) achieved was 0.5 ng/mL, with a linear range of 0.5 ng/mL to 200 ng/mL. Sample recoveries were all greater than 80 % using the external standard method.

**Aeris™ Ordering Information**

PEPTIDE 2.6 μm Minibore Columns (mm)					SecurityGuard™ ULTRA Cartridges*
Phase	50 x 2.1	100 x 2.1	150 x 2.1	250 x 2.1	3/pk
XB-C18	<a href="#">00B-4505-AN</a>	<a href="#">00D-4505-AN</a>	<a href="#">00F-4505-AN</a>	<a href="#">00G-4505-AN</a>	<a href="#">AJ0-8948</a>

for 2.1 mm ID

\*SecurityGuard ULTRA Cartridges require holder, Part No.: [AJ0-9000](#)



**Need a different column size or sample preparation format?**

No problem! We have a majority of our available dimensions up on [www.phenomenex.com](http://www.phenomenex.com), but if you can't find what you need right away, our super helpful Technical Specialists can guide you to the solution via our online chat portal [www.phenomenex.com/Chat](http://www.phenomenex.com/Chat).

**Australia**

t: +61 (0)2-9428-6444  
auinfo@phenomenex.com

**Austria**

t: +43 (0)1-319-1301  
anfrage@phenomenex.com

**Belgium**

t: +32 (0)2 503 4015 (French)  
t: +32 (0)2 511 8666 (Dutch)  
beinfo@phenomenex.com

**Canada**

t: +1 (800) 543-3681  
info@phenomenex.com

**China**

t: +86 400-606-8099  
cninfo@phenomenex.com

**Czech Republic**

t: +420 272 017 077  
cz-info@phenomenex.com

**Denmark**

t: +45 4824 8048  
nordicinfo@phenomenex.com

**Finland**

t: +358 (0)9 4789 0063  
nordicinfo@phenomenex.com

**France**

t: +33 (0)1 30 09 21 10  
franceinfo@phenomenex.com

**Germany**

t: +49 (0)6021-58830-0  
anfrage@phenomenex.com

**Hong Kong**

t: +852 6012 8162  
hkinfo@phenomenex.com

**India**

t: +91 (0)40-3012 2400  
indiainfo@phenomenex.com

**Indonesia**

t: +62 21 3952 5747  
indoinfo@phenomenex.com

**Ireland**

t: +353 (0)1 247 5405  
eireinfo@phenomenex.com

**Italy**

t: +39 051 6327511  
italiainfo@phenomenex.com

**Japan**

t: +81 (0) 120-149-262  
jpinfo@phenomenex.com

**Luxembourg**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**Mexico**

t: 01-800-844-5226  
tecnicomx@phenomenex.com

**The Netherlands**

t: +31 (0)30-2418700  
nlinfo@phenomenex.com

**New Zealand**

t: +64 (0)9-4780951  
nzinfo@phenomenex.com

**Norway**

t: +47 810 02 005  
nordicinfo@phenomenex.com

**Poland**

t: +48 22 51 02 180  
pl-info@phenomenex.com

**Portugal**

t: +351 221 450 488  
ptinfo@phenomenex.com

**Singapore**

t: 800-852-3944  
sginfo@phenomenex.com

**Slovakia**

t: +420 272 017 077  
sk-info@phenomenex.com

**Spain**

t: +34 91-413-8613  
espinfo@phenomenex.com

**Sweden**

t: +46 (0)8 611 6950  
nordicinfo@phenomenex.com

**Switzerland**

t: +41 (0)61 692 20 20  
swissinfo@phenomenex.com

**Taiwan**

t: +886 (0) 0801-49-1246  
twinfo@phenomenex.com

**Thailand**

t: +66 (0) 2 566 0287  
thaiinfo@phenomenex.com

**United Kingdom**

t: +44 (0)1625-501367  
ukinfo@phenomenex.com

**USA**

t: +1 (310) 212-0555  
info@phenomenex.com

🌐 **All other countries/regions  
Corporate Office USA**

t: +1 (310) 212-0555  
[www.phenomenex.com/chat](http://www.phenomenex.com/chat)

**www.phenomenex.com**

Phenomenex products are available worldwide. For the distributor in your country/region, contact Phenomenex USA, International Department at [international@phenomenex.com](mailto:international@phenomenex.com)

**BE-HAPPY™  
GUARANTEE**

Your happiness is our mission. Take 45 days to try our products. If you are not happy, we'll make it right.

[www.phenomenex.com/behappy](http://www.phenomenex.com/behappy)

**Terms and Conditions**

Subject to Phenomenex Standard Terms and Conditions, which may be viewed at [www.phenomenex.com/TermsAndConditions](http://www.phenomenex.com/TermsAndConditions).

**Trademarks**

Aeris and BE-HAPPY are trademarks of Phenomenex. Shimadzu is a registered trademark of Shimadzu Corporation. SCIEX is a registered trademark and Triple Quad is a trademark of AB SCIEX Pte. Ltd.

**Disclaimer**

Phenomenex is in no way affiliated with Shimadzu Corporation.

FOR RESEARCH USE ONLY. Not for use in clinical diagnostic procedures.

© 2024 Phenomenex, Inc. All rights reserved.

